

100 mL Digestion Vessels & Accessories



Our MP-125 digestion tubes are certified to class A volume accuracy, include a metals background report and are compatible with leading brand graphite block digestion systems.

- 250 vessels/case
- Made from virgin polypropylene (PP)
- Excellent visibility
- Very good resistance to acids
- Molded-in graduations at 10, 25, 50, 75 and 100 mL
- Each shipment comes with a certificate of volume accuracy stating the tubes meet Class A tolerance at all graduation lines
- Total capacity of 130 mL
- Extremely low leachable metal content: each shipment comes with a metals background report
- Fit leading brand graphite digestion block systems
- Can be used as a digestion vessel, graduated cylinder, volumetric flask, autosampler tube, or storage container
- Vessels available with or without caps
- Ribbed watch glass available: our item number MP-130

Our 100 mL vessel products fit the following digestion blocks systems:

- EE HotBlock
- SCP Science Digi Block
- LabTech DigiBlock
- Questron Qblock
- Seal Analytical/Thomas Cain SmartBlock II
- Seal Analytical/Thomas Cain Deena II